

**Rapid determination of <sup>90</sup>Sr and <sup>90</sup>Y in water samples with Hidex 300 SL TDCR liquid scintillation counter by Cerenkov-TDCR method**


Introduction

Sr-90 is one of the most hazardous isotopes because of its long physical and biological half-life and because it may cause damage to bone marrow. Therefore rapid and reliable methods to measure radiostrontium are needed particularly in a case of nuclear accident.

Materials and methods

Both Eichrom and Triskem offer excellent products for radiochemical separation of Strontium-90 from its daughter Yttrium-90. Hidex 300 SL automatic TDCR liquid scintillation counter offers excellent measurement capabilities for rapid, robust and easy determination of Strontium quantities. Figure 1 shows flow chart of sample preparation and measurements.

The unique feature of Hidex 300 SL is the TDCR counting which can also be utilized in Cerenkov counting mode.



**Fig. 1** Diagram showing an overview of the methodology used to determine <sup>90</sup>Sr and <sup>90</sup>Y samples.



## Application note

DOC 413-010

HIDEX 300 SL™ TDCR LSC

Version 1

### Procedure

1. Prepare samples by radiochemical purification according to manufacturers instructions. ( Eichrom or Triskem)
2. After radiochemical purification take 8ml sample and measure in 20 ml plastic vial without scintillation cocktail. Count activity of Y-90 directly with Hidex 300 SL TDCR-Cerenkov method. Use readymade parameter file.  
Efficiency Y-90 =  $0.6886 * TDCR + 0.1678$   
\*Use background corrected TDCR value with low activity samples
3. Add 12 ml AquaLight cocktail to sample and count total activity  
Total Bq =  $CPM_{net} / (TDCR * 60)$
4. Subtract results of the first measurement from the results of the second measurement to obtain activities for Sr-90.

### Results and discussion

References customers have reported detection limits of:

#### Sr-90

1 h: 90 mBq/l  
5 h: 40 mBq/l  
10 h: 28 mBq/l  
24 h: 18 mBq/l

Sample volume: 2 L

$k_{1-a} + k_{1-b} = 3.29$

Blank: 50 CPM (for standard model instrument).

Please note that the detection limits are strongly dependent on sample preparation methods.

## Product Information

Finland and International Hidex Oy, Tel. +358 10 843 5570, [info@hidex.com](mailto:info@hidex.com),  
[www.hidex.com](http://www.hidex.com)