



for all your shielding projects.

www.nuclear-shields.com

Hot Lab Cabinets

All configurations



**NUCLEAR
SHIELDS**

Hot Lab Cabinets

Lead lined hot lab cabinets to store radioactive materials. This modular range of shielded hotlab furniture is available with 12.5 or 25 mm lead shielding. Due to the modular design, they can be easily placed next to each other without large open spaces. The only additional thing that has to be ordered is the new size tabletop finishing to prevent leaks. These shielded cabinets are available in multiple configurations, such as three or six drawers. The frames of our shielded hotlab furniture range are made from stainless steel and the shielded front plates are finished with HPL (high pressure laminate) panels, which are often used in clean rooms for their easy decontamination and cleanability.

All cabinets have outside dimensions of 730 x 555 x 996 mm (w x d x h).

Product code	Description	Lead thickness		Inside dimensions (w x d x h)
SC121	Door + retractable shelves (optional)	12.5 / 25	mm	337 x 396 x 152 mm (tray size)
SC123	3 drawers	12.5 / 25	mm	362 x 450 x 116 mm
SC126	6 drawers	12.5 / 25	mm	129 x 450 x 118 mm

SC121

including retractable shelves option



SC126

