



**NUCLEAR  
SHIELDS**

# Syringe shields

size chart



<b>syringe (ml / cc)</b>	<b>Inner diameter (mm)</b>	<b>Inner length (mm)</b>	<b>Big-screw version product code</b>	<b>Clip-lock version product code</b>
1	8.8	72	SS2BS1-1	SS2CL1-1
1	8.8	76	SS2BS1-2	SS2CL1-2
1	11.5	76	SS2BS1-3	SS2CL1-3
2	11.5	53	SS2BS2-4	SS2CL2-4
2	11.5	60	SS2BS2-5	SS2CL2-5
3	11.5	65	SS2BS3-6	SS2CL3-6
3	11.5	69	SS2BS3-7	SS2CL3-7
5	15.4	61	SS2BS5-8	SS2CL5-8
5	15.4	65	SS2BS5-9	SS2CL5-9
5	15.4	68	SS2BS5-10	SS2CL5-10
10	18.2	75	SS2BS10-11	SS2CL10-11
10	18.2	85	SS2BS10-12	SS2CL10-12