

Genie™ 4.0 and Apex® Products

Operating System and Database Qualifications



MIRION
TECHNOLOGIES

Status as of June 2023


Base Model / Name	Version	Legacy Genie Model(s)	Windows 10	Windows 11	Windows Server 2019	Windows Server 2022	SQL Server 2019	SQL Server 2022	Oracle 19c
Genie / Genie™ Spectroscopy Suite	4.0	S500/S502/S504 Genie Basic S501 Gamma Analysis S505 Genie Quality Assurance S506 Genie Interactive Peak Fit S509 Alpha Analysis S561 Genie Batch Programming Genie-FieldPro™ Software	✓	✓	Coming soon	Coming soon			
GenieLite / GenieLite Basic Spectroscopy	4.0	S500/S502/S504 Genie Basic	✓	✓					
ISOCS / ISOCS™/LabSOCS™ Mathematical Efficiency Software	4.5	S573 ISOCS S574 LabSOCS	✓	✓					
Apex-InVivo / Apex InVivo™ Workstation	1.3	S734 Apex-InVivo Workstation	Coming soon	✓			Coming soon		
AIV-Server / Apex-InVivo Server Software	1.3	S738 Apex-InVivo Server				Coming soon	Coming soon		
AIV-Review / Apex-InVivo Review Client	1.3	S735 Apex-InVivo Review Client	Coming soon	Coming soon					
Apex-Alpha / Apex-Alpha™ Desktop/Server	1.4	S770 Apex-Alpha Desktop S772 Apex-Alpha Server	✓		✓		✓		
AA-Client / Apex-Alpha Client	1.4	S771 Apex-Alpha Client	✓						
MGAU / Multi-Group Analysis for Uranium (MGAU™) Software	4.3	S507 MGAU		✓					
MGA / Multi-Group Analysis (MGA™) Software	10.1	S508 MGA		✓					
Apex-Gamma / Apex-Gamma™ Desktop/Server Software	1.5.1	S700 Apex-Gamma Desktop S702 Apex-Gamma Server S703 Apex-Gamma NPP Option	Coming soon	Coming soon	Coming soon	Coming soon	Coming soon	Coming soon	Coming soon
AG-Client / Apex-Gamma Client	1.5.1	S701 Apex-Gamma Client	Coming soon	Coming soon					
Apex-Guard / Apex-Guard™ Option	1.1	S707 Apex-Guard	Coming soon	Coming soon	Coming soon	Coming soon	Coming soon	Coming soon	Coming soon

Product qualification, at time of qualification, for each indicated operating system and configuration is in compliance with the following standards: ISO 9001, ASME NQA-1c-1992, ASME NQA-2-1990 Part 2.7, and latest revision of IEEE 730.1.

Please note: NDA 2000™ software will continue to require Genie 2000 Version 3.4.1 until further notice.

Supported MCAs in Genie Software and Apex Products

The below graphic indicates the years when Mirion MCAs were manufactured (solid bars), and continued support within Genie software after end of MCA manufacturing (thin line).



Genie 4.0 is qualified with the following devices:

- Osprey® DTB
- InSpector™ 2000
- Falcon 5000® Instrument
- DSA-LX® DSA

- Lynx® MCA
- Lynx® II DSA
- Multiport II™ MCA
- Alpha Analyst™ Spectrometer

- Aegis™ Portable HPGe Spectrometer
- Sample Changers

